



SANODURE® VIOLET CLW

Sanodure Violet CLW is quick-striking violet dye suitable for many novelty applications such as costume jewelry, identification tags, decorative appliqués etc.

Chemical description:	Azo dye (granular)
Light fastness:	6 – 7 (good)
Heat fastness:	Moderate up to 300° F.

RECOMMENDED APPLICATION DETAILS:

Concentration:	0.5 - 5.0 g/l
pH range:	5.0 - 6.0
Buffering:	Anodal® Dye Buffer Liquid
Temperature:	140° - 160° F.
Immersion time:	2-20 min.

PREFERRED SEALING METHODS:

Hot:	Anodal ASL or Anodal AS
Medium:	Anodal MS-1 New or Anodal MTS Plus
Cold:	Anodal CS-2, Anodal CS-2N or Anodal CS3

PRECAUTIONS:

The color value of the dye solution is affected by contamination with aluminum and other metal salts. Effective rinsing and care to minimize drag-in are recommended for maximum tank life.

ENVIRONMENTAL CONSIDERATIONS:

This RoHS compliant product contains copper as an integral part of the dye. The preferred method of disposal is described in the Reliant bulletin, the *Anodal WT System*.

Recommendations, notices or instructions as to handling, use, storage or disposal of this product, including its use alone or in combination with other products, or as to any apparatus or process for its use are based upon information believed to be reliable. No liability is taken with respect to any such recommendations or instructions. Sole and exclusive warranty is that products comply with published chemical and physical specifications as provided on the certificate of analysis. No other warranties, either express or implied are given.